

910-3013 Solution 13 – Acridine Orange •DAPI 1 mL

Contents	<p>Microtube containing 1 mL Solution 13 - Acridine Orange (AO) •DAPI</p> <p>The solution contains AO (30 µg/mL ~ 81 µM ~ 0.003 %), 4',6-diamidino-2-phenylindole (DAPI) (100 µg/mL ~ 286 µM ~ 0.01 %) and sodium azide (NaN₃) (0.1 mg/mL ~ 1.5 mM ~ 0.01 %).</p> <p>EINECS no. for AO is 233-353-6 EINECS no. for DAPI is 249-186-7 EINECS no. for NaN₃ is 247-852-1</p>
Application	<p>Solution 13 - AO •DAPI is useful for measuring cell concentrations and viability of nucleated animal cells.</p>
Principle	<p>Solution 13 - AO •DAPI contains two different dyes: AO staining all nucleated cells and DAPI staining nonviable cells only.</p>
Usage	<p>The solution is for research and development purposes only and is not for diagnostic or therapeutic use.</p>
Storage	<p>Store at 2-7°C</p>
Stability	<p>For unopened microtubes, the expiry date is shown on the microtube and on the plastic container. The solution is produced 15 months before the expiry date. The solution expires 2 months after opening the microtube.</p>
Safety Information	<p>In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. If a spill is observed, perform a clean-up of the area, which may have been in contact with the solution. Use gloves and suitable protective clothing. Please also refer to MSDS regarding safety information.</p>
Disposal of Waste	<p>After use, the microtube should be disposed of according to national or local laws and regulations regarding the nature of the mixture it contains.</p>